
GB

LIFECATH-TWIN TWO SEPARATE POLYURETHANE CUFFED DIALYSIS CATHETERS INSTRUCTIONS FOR USE Code: 1276., 8276.

PRODUCT DESCRIPTION The LIFECATH-TWIN is a radiopaque polyurethane catheter with a lockable extension and a Dacron cuff.

INDICATIONS The polyurethane LIFECATH-TWIN is designed for acute and chronic haemodialysis. LIFECATH-TWIN is percutaneously inserted via the internal jugular, subclavian or femoral veins in the adult patient.

CONTRAINDICATIONS Do not use this catheter in thrombosed veins.

- POTENTIAL COMPLICATIONS - Sepsis - Exsanguination - Central vein thrombosis - Subcutaneous hematoma - Trauma to major vessel or right atrium - Thrombosis - Pulmonary emboli - Haemorrhage - Haemothorax - Exit site infection - Vessel trauma - Cardiac tamponade - Cardiac arrhythmia - Subcutaneous tunnel infection - Haematoma - Luminal thrombosis - Pneumothorax

CATHETER PRECAUTIONS & WARNINGS Do not restyle the catheter or accessories by any method.

Contents sealed by re-use or resterilisation of the catheter or accessories. Contaminated and non-sterile catheter in unopened or undamaged package. Do not use catheter if package is damaged or has been opened. Do not use catheter if components show signs of damage (crimped, crushed, etc).

Before manipulating the insertion, ensure that you are familiar with the above complications and their emergency treatment should any of them occur. In the rare event that a hub or connector separates from any component during the insertion or use, take all necessary steps and precautions to prevent blood loss or air embolism and remove catheter immediately.

Do not advance the stainless steel guidewire or catheter if unusual resistance is encountered. Do not insert or withdraw the guidewire forcibly from any component. The guide wire could break or unravel, in which case both the catheter and guidewire must be removed simultaneously.

Do not use sharp instruments near the extension lines or tubing. Do not use scissors to remove dressing, as this could possibly cut or damage the catheter. Do not suture through any part of the catheter. Catheter tubing can tear when subjected to excessive force or rough edges.

Use only smooth jawed forceps for clamping when using the clamp supplied with the catheter. We recommend using only fine extension clamps which have been provided for clamping. Clamping the catheter repeatedly in the same spot could weaken the tubing.

It is recommended that only laser lock connections be used with this catheter (including syringes, bloodlines, IV tubing, and injection caps). Repeated over tightening of bloodlines, syringes, and caps will reduce contact area and could lead to potential connector failure.

Catheter implantation duration shall be determined by the physician based on catheter patency, infection, and indications for use. The medical techniques and procedures described in these instructions do not represent all medically acceptable protocols.

Puncture of the right atrium or superior vena cava with vessel dilator or introducer may result in serious injury or death. Prolonged exposure to ultraviolet light may damage the catheter.

CATHETER INSERTION DIRECTIONS INSERTION SITE Although the polyurethane catheter may be inserted in the subclavian, jugular or femoral, the right internal jugular is the preferred insertion site.

And while the subclavian vein is superior for patient comfort and ease in changing dressings, long term use of this vein with thrombus formation and vein stenosis.

PREPARATION The operation theatre is the preferred location for the procedure; however bedside insertion is acceptable if transport to the operating theatre is not possible and if the procedure is performed under sterile conditions.

1. Provide a sterile operative field; wear sterile drapes, instrument and accessories. Perform surgical scrub. Wear gown cap, gloves and mask.

2. Shave the access site (optional) and prep the site with an appropriate iodine-free antiseptic and wait to dry. Isolate the access site with sterile drapes.

3. Administer a local anesthetic to the skin and underlying tissue.

4. Prime the catheter with sterile heparinised saline and clamp the extension tubes. WARNING: To prevent air embolism keep the catheter clamped at all times when not attached to a syringe or bloodline.

PERCUTANEOUS INSERTION METHOD CREATING THE SUBCUTANEOUS TUNNEL SELDINGER INSERTION METHOD Anesthetise the insertion site, attach a syringe to the introducer needle. Insert needle into identified vein of choice and aspirate to ensure correct placement.

5. Prime the catheter with sterile heparinised saline and clamp the extension tubes. WARNING: To prevent air embolism keep the catheter clamped at all times when not attached to a syringe or bloodline.

6. CAUTION: Monitor the patient during the procedure. Cardiac arrhythmias may result if the guidewire is allowed to pass into the right atrium. Hold the guidewire in place and remove the needle (Fig. 3).

7. CAUTION: Never pull the guidewire back through the needle, this could result in serious damage to the guidewire. Advance the split sheath introducer over the guidewire, dilate the insertion tract (at the point of skin entry, you may wish to widen the insertion tract by the use of a scalpel).

8. Remove the dilator/guidewire from the split sheath introducer. Prime the LIFECATH-TWIN catheter, insert the catheter through the sheath and into the vein. Hold the catheter firmly in place and withdraw the sheath from the vein by pulling both tabs of the sheath out of the entry site simultaneously at a 180 degree angle.

9. The sheath will separate from the catheter, leaving the catheter in the vein. Verify catheter placement by x-ray.

10. Place the temporary colour coded clamping sleeves with clamps onto each catheter to prevent bleed back or the risk of air embolism.

11. Position catheter over desired tunnel length. Identify the desired point at which the cuff will be positioned to pass into the right atrium.

12. Anesthetise the entire length of the tunnel. Make an incision at the tunnel exit site. Connect the trocar and football tunneller together using the tunnelling devices make an 8-10cm tunnel in the direction of the catheter exit site.

13. The football tunneller should only dilate up to the point where the cuff of the catheter is intended to sit. To help widen this part of the tunnel slide the tunneller up and down the cuff about 6-8 times. Remember the wider area should only be up to the point where the cuff will be positioned. This will create a stop for the cuff so it cannot migrate out of the insertion site or exit site.

14. Pull the football tunneller out of the tunnel through the catheter insertion site. Connect the proximal end of the catheter onto the trocar and remove the clamp on the catheter.

15. Carefully pull the catheter through the tunnel feeding the cuff into the tunnel as well. Palpate the tunnel until the cuff is situated in the desired position. The catheter may be cut to the desired length and the Lifecath Lock hubs attached.

Repeat the 11) to 15) for the second catheter.

TIP: A tunnel with an arc reduces the possibility of the catheter kinking in the tract and resulting in poor blood flow.

NOTE: If arterial blood is aspirated, remove the needle and apply immediate pressure to the site for at least 15 minutes. Ensure that the bleeding has stopped and that no hematoma has developed before attempting to cannulate the vein again.

CAUTION: Failure to insert the wire promptly may lead to blood loss through the needle. Cardiac arrhythmia may result if the guidewire passes into the right atrium. If symptoms occur, pull back the guidewire.

WARNING: Puncture of the vessel or right atrium with the vessel dilator may result in serious injury or death. Ensure that the guidewire does not move further into the vein. Hold the sheath in place and gently remove the dilator and guidewire and discard both.

WARNING: To prevent air embolism and excessive bleeding, pinch the sheath closed as the dilator is removed. CAUTION: Do not allow the catheter to move out of the vein with the sheath. Ensure that the vein is not bleeding around the catheter. Use fluoroscopy or portable x-ray to view the catheter. Position the catheter in the superior vena cava.

WARNING: To prevent air embolism expel all air from syringe prior to injecting saline and heparin. Inject 10 ml of sterile normal saline into each lumen, then quickly inject heparinised saline using positive flush technique. Repeat until the catheter is free of air. Repeat this procedure until the catheter is free of air.

CAUTION: Remove skin suture by the 4th or 5th day to diminish erythema in the area. The catheter can also be immobilised with either a bandage and gauze or transparent dressings.

Attach Lifecath Lock, see Fig 4 to 10. Fig 4: Illustrates the complete compression device. Fig 5: Remove the white protective plastic sleeve. Fig 6: Insert the metal rod into the catheter (ensuring green Lifecath Lock sleeve is still in situ) then slide Lifecath Lock hub onto catheter.

Fig 7: Slide catheter rod until the green Lifecath Lock is covered and the line is completely flush with the main Lifecath Lock body. Slide the Lifecath Lock hub over the catheter to meet the main body. Fig 8: Push the securement hole with the securement locator and push until an audible click is heard. Fig 9: Remove metal rod ensuring line is clamped post removal.

Fig 10: Attach a syringe to the extension and open the clamps, this should allow blood to aspirate. If any resistance to aspiration is observed, the catheter may need to be rotated or repositioned, to allow good blood flow. Flush the extension with heparinised saline in accordance with your Hospital policy, replace the clamp and attach the injection membrane. In order to avoid contamination, such as in the rare event of a leak, the catheter should be clamped to prevent blood loss and the physician/health care professional consulted.

CAUTION: Clamp only on the extension tubing, never on the actual catheter. Do not use serrated forceps, use only the in-line clamps provided. Confirm placement of the catheter immediately after insertion by use of X-ray. Observe the patient for any signs of cardiac arrhythmia. Record inflowing catheter length.

CAUTION: Failure to verify catheter placement may result in serious trauma or fatal complications. Frequent visual inspections should be carried out before delays begin to ensure the integrity of the catheter and dilator. In the rare event of a leak, the catheter should be clamped to prevent blood loss and the physician/health care professional consulted.

When removing the catheter, DO NOT use a sharp, jerking motion or undue force, this may damage the catheter. The cuff should be freed from the tissue prior to removal. After removal manual pressure should be applied to the puncture site for 10-15 minutes until there is no further bleeding. A transparent wound dressing should be applied.

WARNING: If the catheter offers resistance do not pull further. Perform a cutdown removal all sutures at the venotomy site.

HEPARINISATION 1. Prime the catheter with sufficient volume and concentration of heparin to maintain patency. Use the heparin concentration approved by your hospital policy.

2. Aspirate inflowing heparin from each catheter before infusing fresh heparin. After aspirating, flush each catheter thoroughly with sterile, normal saline before infusing fresh heparin. Attach sterile injection caps to the luer connectors.

3. Clamp immediately after infusing to create positive pressure in the lumens. The techniques and procedures described in these instructions do not represent all medically acceptable protocols.

4. Once the catheter has been primed with heparinised saline, do not open the clamp of the extension tube. If the clamp is opened, the lumene volume will increase slightly as a result of the extension tube sucking into its "normal" unclamped state. This creates a vacuum at the tip, causing blood to be returned to the distal portion of the catheter, ultimately resulting in a thrombus formation (positive central venous pressure).

5. If the catheter is used immediately after insertion, the concentration of heparin-saline used to flush the catheter can be relatively low.

MANAGEMENT OF ONE-WAY OBSTRUCTION One-way obstruction exists when a lumen can be flushed but cannot be aspirated. Reversing the bloodlines may resolve the obstruction using this method as an alternative to catheter replacement. Dialyse the patient by connecting the arterial bloodline to the venous adaptor and the venous bloodline to the arterial adaptor. A significant increase in recirculation may be expected.

THROMBUS FORMATION Never forcibly flush an obstructed lumen. If either lumen develops a thrombus, first attempt to aspirate the clot with a syringe. If aspiration fails the physician/health care professional may attempt to dissolve the clot with urokinase.

WARNING: Urokinase may cause systemic fibrinolysis if infused into the circulation. Refer to manufacturers instructions, indications for use and contraindications before using urokinase. Streptokinase is not recommended; it is reported to be highly anaphylogenic.

DIALYSIS Hepatitis analysis machine and prime bloodlines in the usual manner. Maintain sterile technique when handling or using the catheter. Scrub the adaptors, injection caps, clamps and extension tubes of the catheter with a disinfection solution.

2. Verify that the arterial extension tube is clamped before removing the injection cap. Remove the injection cap from the adaptor and attach a luer lock syringe. Confirm patency of the lumen by aspirating inflowing heparinised saline and then aspirating the lumen. Heparinised saline is used if may be administered at this time. Leave the syringe in place. Clamp and close the extension tube.

3. Repeat step 2 for the venous lumen.

WARNING: When connecting bloodlines to the catheter do not permit air to enter the blood path. Remove the syringe and connect the arterial bloodline to the arterial adaptor, open the clamp on the arterial extension and turn on the blood pump.

5. Prime the extracorporeal circuit with the patient's blood and turn off the blood pump. Ensure the venous extension tube is clamped, then remove the syringe and connect the venous bloodline to the venous adaptor. Open the clamp on the venous extension tube and turn the pump on.

POST DIALYSIS 1. Observe sterile technique. Stop the blood pump. Close the clamp on the arterial bloodline at the connection site. Disconnect the arterial bloodline from the adaptor of the catheter.

2. Connect a 10 ml syringe filled with normal saline to the arterial adaptor, open the clamp on the arterial extension tube and flush the blood from the arterial lumen of the catheter.

3. Connect saline solution to the venous adaptor, open the saline clamp and flush back the patient's blood through the venous lumen of the catheter.

4. After the patient's blood has been flushed back, turn off the blood pump. Clamp the venous extension tube and disconnect the venous bloodline to the venous adaptor of the catheter. Connect a 10 ml syringe filled with appropriate concentration of heparin solution to the venous adaptor, open the clamp on the venous extension tube and quickly infuse the heparin. Reclamp the extension tube, removing the syringe and attach a sterile injection cap to the adaptor.

5. Verify that the arterial extension tube is clamped. Connect a 10 ml syringe containing the appropriate concentration of heparin solution to the arterial adaptor of the catheter, open the clamp on the arterial extension tube and quickly infuse the heparin. Reclamp the extension tube, remove the syringe and attach a sterile injection cap to the adaptor.

WARNING: Keep the catheter clamped at all times except when connected to the blood lines or syringe during treatment. EXIT SITE CARE Exit site care is required at each treatment. Use sterile technique when handling or using the catheter. NOTE: Wound dressings must be kept dry. The patient must not swim, shower or wet the dressing while walking. If adhesion is compromised by perspiration or accidental wetting, the dressing must be changed by the qualified medical personnel under sterile conditions.

CAUTION: Failure to insert the wire promptly may lead to blood loss through the needle. Cardiac arrhythmia may result if the guidewire passes into the right atrium. If symptoms occur, pull back the guidewire.

WARNING: Puncture of the vessel or right atrium with the vessel dilator may result in serious injury or death. Ensure that the guidewire does not move further into the vein. Hold the sheath in place and gently remove the dilator and guidewire and discard both.

WARNING: To prevent air embolism and excessive bleeding, pinch the sheath closed as the dilator is removed. CAUTION: Do not allow the catheter to move out of the vein with the sheath. Ensure that the vein is not bleeding around the catheter. Use fluoroscopy or portable x-ray to view the catheter. Position the catheter in the superior vena cava.

WARNING: To prevent air embolism expel all air from syringe prior to injecting saline and heparin. Inject 10 ml of sterile normal saline into each lumen, then quickly inject heparinised saline using positive flush technique. Repeat until the catheter is free of air. Repeat this procedure until the catheter is free of air.

CAUTION: Remove skin suture by the 4th or 5th day to diminish erythema in the area. The catheter can also be immobilised with either a bandage and gauze or transparent dressings.

Attach Lifecath Lock, see Fig 4 to 10. Fig 4: Illustrates the complete compression device. Fig 5: Remove the white protective plastic sleeve. Fig 6: Insert the metal rod into the catheter (ensuring green Lifecath Lock sleeve is still in situ) then slide Lifecath Lock hub onto catheter.

Fig 7: Slide catheter rod until the green Lifecath Lock is covered and the line is completely flush with the main Lifecath Lock body. Slide the Lifecath Lock hub over the catheter to meet the main body. Fig 8: Push the securement hole with the securement locator and push until an audible click is heard. Fig 9: Remove metal rod ensuring line is clamped post removal.

Fig 10: Attach a syringe to the extension and open the clamps, this should allow blood to aspirate. If any resistance to aspiration is observed, the catheter may need to be rotated or repositioned, to allow good blood flow. Flush the extension with heparinised saline in accordance with your Hospital policy, replace the clamp and attach the injection membrane. In order to avoid contamination, such as in the rare event of a leak, the catheter should be clamped to prevent blood loss and the physician/health care professional consulted.

CAUTION: Clamp only on the extension tubing, never on the actual catheter. Do not use serrated forceps, use only the in-line clamps provided. Confirm placement of the catheter immediately after insertion by use of X-ray. Observe the patient for any signs of cardiac arrhythmia. Record inflowing catheter length.

CAUTION: Failure to verify catheter placement may result in serious trauma or fatal complications. Frequent visual inspections should be carried out before delays begin to ensure the integrity of the catheter and dilator. In the rare event of a leak, the catheter should be clamped to prevent blood loss and the physician/health care professional consulted.

When removing the catheter, DO NOT use a sharp, jerking motion or undue force, this may damage the catheter. The cuff should be freed from the tissue prior to removal. After removal manual pressure should be applied to the puncture site for 10-15 minutes until there is no further bleeding. A transparent wound dressing should be applied.

WARNING: If the catheter offers resistance do not pull further. Perform a cutdown removal all sutures at the venotomy site.

HEPARINISATION 1. Prime the catheter with sufficient volume and concentration of heparin to maintain patency. Use the heparin concentration approved by your hospital policy.

2. Aspirate inflowing heparin from each catheter before infusing fresh heparin. After aspirating, flush each catheter thoroughly with sterile, normal saline before infusing fresh heparin. Attach sterile injection caps to the luer connectors.

3. Clamp immediately after infusing to create positive pressure in the lumens. The techniques and procedures described in these instructions do not represent all medically acceptable protocols.

4. Once the catheter has been primed with heparinised saline, do not open the clamp of the extension tube. If the clamp is opened, the lumene volume will increase slightly as a result of the extension tube sucking into its "normal" unclamped state. This creates a vacuum at the tip, causing blood to be returned to the distal portion of the catheter, ultimately resulting in a thrombus formation (positive central venous pressure).

5. If the catheter is used immediately after insertion, the concentration of heparin-saline used to flush the catheter can be relatively low.

MANAGEMENT OF ONE-WAY OBSTRUCTION One-way obstruction exists when a lumen can be flushed but cannot be aspirated. Reversing the bloodlines may resolve the obstruction using this method as an alternative to catheter replacement. Dialyse the patient by connecting the arterial bloodline to the venous adaptor and the venous bloodline to the arterial adaptor. A significant increase in recirculation may be expected.

THROMBUS FORMATION Never forcibly flush an obstructed lumen. If either lumen develops a thrombus, first attempt to aspirate the clot with a syringe. If aspiration fails the physician/health care professional may attempt to dissolve the clot with urokinase.

WARNING: Urokinase may cause systemic fibrinolysis if infused into the circulation. Refer to manufacturers instructions, indications for use and contraindications before using urokinase. Streptokinase is not recommended; it is reported to be highly anaphylogenic.

DIALYSIS Hepatitis analysis machine and prime bloodlines in the usual manner. Maintain sterile technique when handling or using the catheter. Scrub the adaptors, injection caps, clamps and extension tubes of the catheter with a disinfection solution.

2. Verify that the arterial extension tube is clamped before removing the injection cap. Remove the injection cap from the adaptor and attach a luer lock syringe. Confirm patency of the lumen by aspirating inflowing heparinised saline and then aspirating the lumen. Heparinised saline is used if may be administered at this time. Leave the syringe in place. Clamp and close the extension tube.

3. Repeat step 2 for the venous lumen.

WARNING: When connecting bloodlines to the catheter do not permit air to enter the blood path. Remove the syringe and connect the arterial bloodline to the arterial adaptor, open the clamp on the arterial extension and turn on the blood pump.

5. Prime the extracorporeal circuit with the patient's blood and turn off the blood pump. Ensure the venous extension tube is clamped, then remove the syringe and connect the venous bloodline to the venous adaptor. Open the clamp on the venous extension tube and turn the pump on.

POST DIALYSIS 1. Observe sterile technique. Stop the blood pump. Close the clamp on the arterial bloodline at the connection site. Disconnect the arterial bloodline from the adaptor of the catheter.

2. Connect a 10 ml syringe filled with normal saline to the arterial adaptor, open the clamp on the arterial extension tube and flush the blood from the arterial lumen of the catheter.

3. Connect saline solution to the venous adaptor, open the saline clamp and flush back the patient's blood through the venous lumen of the catheter.

4. After the patient's blood has been flushed back, turn off the blood pump. Clamp the venous extension tube and disconnect the venous bloodline to the venous adaptor of the catheter. Connect a 10 ml syringe filled with appropriate concentration of heparin solution to the venous adaptor, open the clamp on the venous extension tube and quickly infuse the heparin. Reclamp the extension tube, removing the syringe and attach a sterile injection cap to the adaptor.

5. Verify that the arterial extension tube is clamped. Connect a 10 ml syringe containing the appropriate concentration of heparin solution to the arterial adaptor of the catheter, open the clamp on the arterial extension tube and quickly infuse the heparin. Reclamp the extension tube, remove the syringe and attach a sterile injection cap to the adaptor.

WARNING: Keep the catheter clamped at all times except when connected to the blood lines or syringe during treatment. EXIT SITE CARE Exit site care is required at each treatment. Use sterile technique when handling or using the catheter. NOTE: Wound dressings must be kept dry. The patient must not swim, shower or wet the dressing while walking. If adhesion is compromised by perspiration or accidental wetting, the dressing must be changed by the qualified medical personnel under sterile conditions.

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