**MultiCath/ MultiStar/ SELDIPUR/ MultiCath Expert/ LeaderCath Expert**  
**Codice 155., 157., 158., 159./ 6155., 6158., 6159./ 128./ 8155., 8157., 8158., 8159., 8128.**

Instructions for use


**Before inserting the catheter:**  
 1. Prep and drape area surrounding intended venipuncture site, create skin wheal at incision site with local anesthetic agent. To prevent air embolism and to facilitate venous puncture, place patient in an appropriate Trendelenburg position.

2. Assess the length of the catheter that needs to be introduced; centimeter markings on the catheter assist in correct placement.


3. Identify the different lumens and corresponding hubs of the catheter. Match intended use of lumens to patient's clinical requirements, e.g. smallest lumen for low delivery rates. Flow rates for individual lumens of different models are shown on the product pack.

4. All lumens other than the guidewire lumen (distal hub) should be flushed, primed with normal or heparinized saline and kept closed with the injection caps or stopcocks provided.  
**Catheter introduction:**  
 When the "J" guidewire is not presented in an Advancer device, the "J" wire with straightener should be removed from the storage coil and used in the conventional manner, taking into account all cautions and warnings.


5. Retract "J" of wire within protector and remove cap.


6. Attach needle to syringe, using a negatively charged syringe. Perform venipuncture. Successful venipuncture is confirmed by free flow of blood into the syringe.

7. Remove syringe from needle and attach advancer, using the thumb advance the guidewire stage by stage until the 20 cm mark is within the needle hub. Gently remove advancer from needle hub, hold guidewire and remove advancer from wire. Remove needle. Adjust position of wire until 20 cm mark can be seen at the chest exit site. Place catheter over wire. Ensure there is sufficient guidewire available so as to hold the wire at the proximal end of the distal hub before proceeding to advance the catheter into the vein.


8. Gently remove advancer from needle hub, hold wire in position and totally remove advancer from wire.


**Caution:**  
 Never pull back or withdraw the guidewire while the needle is in place. This could damage the wire on the needle level. If any resistance is felt to the advancement of the guidewire, do not persist, withdraw both the needle and the guidewire simultaneously.  
 To further decrease the risk of damage to the guidewire, a short i.v. cannula may be used as an introducer.

9. Holding the guidewire in place, remove the needle.


GB USA

10. If appropriate, the venipuncture site may be enlarged with a scalpel. Take care not to cut the guidewire.


11. Dilate the entry point into the vessel by introducing the vessel dilator over the guidewire.


**CAUTION:**  
 When introducing a catheter or dilator over the guidewire make sure that the length of the exposed guidewire is sufficient so that you can maintain a firm grip on the guidewire at all times. The vessel dilator should not be left in the vessel as an indwelling catheter, to do so should cause vessel damage.

12. After removal of dilator, introduce the distal tip of the catheter over the guidewire.


To aid smooth insertion of the catheter hold the catheter close and parallel to the skin and progressively advance catheter to the required depth.

13. Holding the catheter in position at required depth, carefully remove guidewire.




If resistance is felt during guidewire removal, the tip of the guidewire may be kinked within the vessel. Do not persist, but withdraw the catheter relative to the guidewire about 2-3 cm and continue to remove the guidewire. If you still feel resistance withdraw both the catheter and the guidewire simultaneously.

Attach a syringe to main lumen hub and aspirate to check that a free flow of blood confirms correct positioning.


14. Connect main lumen, and any other lumens required at this point to appropriated and prepared equipment.

15. Secure the catheter to the patient using the integral side wings as a primary suture site. The flexible fixation wing should be used as a secondary fixation site if necessary. Spread the wing apart and mount on the catheter from below. The wing position can be changed by again flexing the wing.

Federal U.S. law restricts this device to sale by or on the order of a physician.

16. Attach fastener to the wing over the wing groove and press firmly together.

17. The wing can then be sutured to the skin. Do not suture directly around the catheter as this may compromise the lumens.

18. Check the catheter and dressing regularly for evidence of catheter movement, interrupted flow rates, leakage, security of all connections and signs of infection. Particular care should be taken when changing lines and dressings that the procedure does not damage the catheter. Do not kink or use forceps to clamp the lumen tubings during line changes (use the provided stopcocks instead), avoid potential problem of cutting catheter by not using scissors to remove dressings, etc... The clamps on the lead should only be used temporarily.

19. Nursing information available upon request. Ask for reference DI.08.

**CAUTIONS: Read before inserting and using the catheter!**

1. Correct positioning of the catheter should always be confirmed by x-ray. Correct position will depend on clinical application and introduction point of catheter.

2. The above problems highlight the importance of selecting the most appropriate length of catheter for a particular patient. Left sided neck approaches normally require a longer catheter, e.g. 30 cm length. If longer catheters are used for a right sided neck approach, particular care must be taken not to advance them too far and to ensure there is no accidental migration of the catheter after placement. Always confirm correct placement by x-ray and observe guidelines on assessing required catheter length and the use of catheter distance markings for both placement and checking for migration. The catheter is graduated from its distal tip at 10 mm increments. The first of these is at 90 mm and the last at 160 mm, 200 mm or 300 mm, according to overall catheter length.

3. Reported complications with catheters are extensive and include hydromediastinum, hydrothorax, catheter embolism, air embolism, thrombosis, laceration of the thoracic duct, temporary paralysis of the phrenic nerve, diaphragmatic paralysis, knotting of the catheter, arterial puncture, bacteremia and septicemia. Follow hospital or departmental protocol to prevent air embolism in catheter placement and subsequent maintenance. Observe strict aseptic technique during placement, maintenance and care of catheter.

4. Be aware that organic solvents such as alcohol or acetone may interact with catheter material and weaken it.

5. Only use luer lock connecting tubing with this product.

6. Infusion pressure applied through this catheter should not exceed 1.5 bar.

7. This product may only be sold to or by order of a physician.

Do not use this product for diagnosis, control or treatment of a defect of the heart or the central venous system.

**IMPORTANT CAUTIONS:**  
 The using physician should be familiar with and undertake all standard precautions to avoid problems or complications associated with the introduction and use of sheath introducers. Complications associated with percutaneous sheath introducers include vessel wall perforation, cardiac tamponade, haemorrhage, pleural and mediastinal injury, air embolism, catheter embolism, thoracic duct injury, septicemia, bacteraemia, thrombosis and accidental arterial puncture. The dilator is not suitable for this purpose and must be removed following sheath introduction.

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