




NUTRILINE Long-term catheter with splitting needle GB

Instructions for use

Code 1252, 1262, 1253, 1263, 1254, 1264

The catheter is made of radiopaque PUR, with markings at each 50 mm. It has its own integral extension set, complete with fixation wings. The Luer-Lock connecting adapter at the end of the catheter allows safe connection of an infusion tube. The catheter is introduced through a splitting needle, which comes with a needle hub adaptor and syringe for safe, clean, blood aspiration.

The catheter is a PIC line (Peripherally Inserted Catheter). The catheter is used for mid to long term infusion of TPN or intravenous medications, it is also suitable for peripheral TPN in adults.

Introduction is generally via the basilic, median cubital, external jugular, long saphenous, or, in neonates, the temporal or femoral veins.

Strict aseptic technique must be observed. In order to avoid air embolism the catheter must be flushed before it is introduced.

Caution:
Do not stretch the catheter or subject it to pressure above 21,75psi (1,5 bar, 1,140mmHg). **Do not use syringes smaller than 10 cc, as these can generate very high pressures.** It is possible to generate 4 or 5 times the maximum safety pressure, with any size of hand held syringe. Subjecting the catheter to pressure above 21,75psi can result in catheter rupture and embolism.

Smaller syringes generate higher pressures than larger ones.
Assess the patient and determine the site for introduction of the catheter, and, if cutting the catheter to length, measure to determine the length required. Wash hands thoroughly, don gloves, and rinse gloves free of starch with sterile water.

Carefully clean the insertion site and position drapes. The puncture may only be made after the disinfectant has been completely dried on the skin, or do not clean the tubing with alcohol, acetone or spirit based cleaning solutions as these can damage the tubing. Aqueous based solutions are recommended. Trim catheter to length required, and flush with saline.

Connect the needle hub adaptor into the needle hub, and connect to the syringe. Perform venipuncture using the splitting needle provided.

Do not grip the two wings tightly, as this may cause the needle to break prematurely. Venipuncture is confirmed by a free flowing blood flash back.

Remove the needle hub adaptor and syringe from the introducing needle.

Fig. 1: Puncture with splitting needle.
Insert the Catheter through the needle, using non-toothed forceps. Advance the catheter using short, steady strokes.

Fig. 2: Feeding the catheter with an anatomical forceps.
If difficulty is experienced advancing the catheter, gentle flushing of the catheter with saline solution, may serve to stiffen the catheter slightly, and may also assist in opening venous valves which may be obstructing advancement. When the final position is reached, stabilize the catheter by pressing lightly on the skin, 2-3 inches above the insertion site and withdraw the needle over the catheter.

INTRODUCER REMOVAL.
The introducer is removed from the catheter in five stages:
Stabilize the catheter to prevent withdrawing the catheter with the needle.

IMPORTANT CAUTION: AT NO TIME SHOULD THE CATHETER BE WITHDRAWN BACK THROUGH A SPLITTING NEEDLE. IF IT BECOMES IMPOSSIBLE TO ADVANCE THE CATHETER INTO A SATISFACTORY POSITION, THEN THE NEEDLE AND CATHETER MUST BE WITHDRAWN SIMULTANEOUSLY. THE RESULT OF WITHDRAWING A CATHETER BACK THROUGH THE NEEDLE CAN BE CATHETER EMBOLISM.

Fig. 3: withdraw the needle from the vessel. Do not split or peel away the needle until it is completely removed from the vessel.

Fig. 4: gently pinch the wings of the needle completely together to break the bottom and the top part part of the needle hub.

Fig. 5: peel the wings of the needle apart, until the needle has peeled to within app. 1mm of the tip. The top part of the needle will be completely open, and the catheter can be carefully removed with the forceps. The needle can then be recapped using the one-handed technique, for safe disposal.

Note: If considered preferable, the needle can be peeled completely into two pieces, by increasing, slightly, the peeling force, when resistance is felt, after peeling the needle to 1 mm from the tip

Anchor the catheter, using the fixation wings. If a vein in the arm is used, a small loop should be left in the catheter, to prevent kinking of the catheter if the patient flexes or bends the arm.

IMPORTANT: ARRANGE FOR A CHECK X-RAY TO CONFIRM CORRECT PLACEMENT PRIOR TO STARTING ANY COURSE OF TREATMENT VIA THE CATHETER.

If the catheter is not to be used immediately, or its use is temporarily discontinued, it should be heparin locked to prevent clotting or blockage. CAUTION: Do not over stretch the catheter as it may rupture, and rebound into the insertion site, causing a catheter embolism.

The catheter is not suitable for blood aspiration.

Caution:
Do not exceed a bolus injection pressure of 1,5 bar.
Do not exceed an infusion pressure of 1,5 bar as this will burst the catheter.

Check before starting an infusion, with the help of an X-ray photograph or during placement with the help of an intravascular electrocardiographic lead, that the catheter is in the correct position and record the latter. The catheter tip must not be advanced into the heart (right atrium). The location of the catheter within the heart may cause cardiac tamponade, or cardiac arrhythmias.

Do not use this product for diagnosis, control or treatment of a defect of the heart or the central venous system.

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