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




PREMICATH
Catheter with splitting needle - Code 1261

INSTRUCTIONS FOR USE

The catheter is made of radiopaque PUR, with thin line markings at every 10 mm, two thin line marks at 50 mm and 150 mm, and one thick line mark at 100 mm and 200 mm. A black mark is also found on the distal tip for user to confirm that the entire catheter has been withdrawn (for untrimmed catheter). It has its own integral extension set, complete with fixation wings. The Luer-Lock connecting adaptor at the end of the catheter allows safe connection of an infusion tube. The catheter is introduced through a breakaway needle, peelable cannula or Neocath Split (see separate instructions for use for Neocath Split and peelable cannula enclosed), all of which come with a needle hub adaptor and syringe for safe, clean, blood aspiration.

The catheter is a PICC line (Peripherally Inserted Central Catheter). The catheter is used for mid to long term infusion of TPN or intravenous medications. Introduction is generally via the basilic, median cubital, jugular, long saphenous, temporal or femoral veins.

Strict aseptic technique must be observed. In order to avoid air embolism the catheter must be flushed before it is introduced.

Assess the patient to determine the site for introduction of the catheter. If it is necessary to shorten the catheter, use a scalpel for this purpose to produce a vertical cut, record trimmed catheter length in patient notes. **CAUTION: Prior to trimming the Premicath with Stylet, the stylet must be retracted inside the catheter so that the tip is 1 cm above the anticipated point of trimming. In order to retract stylet, disconnect between the T-Piece and catheter hub, do not try and remove the cap. Do not cut the stylet. Bend the stylet to at least 90° at the catheter hub to prevent stylet from protruding beyond the catheter tip. Any complications resulting therefrom do not lie within the scope of responsibility of the manufacturer.**

Wash hands thoroughly, don gloves, and rinse gloves free of starch with sterile water.

Carefully clean the insertion site and position drapes. The puncture may only be made after the disinfectant has been completely dried on the skin. Be aware that organic solvents such as alcohol or acetone may interact with catheter material and weaken it.

Trim catheter to length if required and flush with saline (**CAUTION: Prior to trimming the Premicath with Stylet, the stylet must be retracted inside the catheter 1 cm from the anticipated point of trimming. In order to retract stylet, disconnect between the T-Piece and catheter hub, do not try and remove the cap. Do not cut the stylet. Bend the stylet to at least 90° at the catheter hub to prevent stylet from protruding beyond the catheter tip.**)

Caution: Do not stretch the catheter or subject it to pressure above 21,75 psi (1,5 bar, 1,140 mmHg). **Do not use small syringes as these can generate very high pressures.** It is possible to generate 4 or 5 times the maximum safety pressure, with any size of hand held syringe. Subjecting the catheter to pressure above 21,75 psi can result in catheter rupture and embolism. Smaller syringes generate higher pressures than larger ones.

Fig. 1 Perform venipuncture using the breakaway needle provided.

NOTE: See individual instructions for use for Neocath Split and Peelable cannula.

Do not grip the two wings tightly, as this may cause the breakaway needle to break prematurely. Venipuncture is confirmed by a free flowing blood flash back.

Caution: Do not apply direct pressure over the catheter if there is a guidewire in situ during insertion. This may cause the guidewire to become kinked and will make it extremely difficult to withdraw the guidewire from the catheter.

Stabilize catheter at the hub and slowly remove guidewire from the catheter using a steady, gentle force.

Caution: In order to remove the stylet, disconnect between the T-Piece and catheter hub, do not try and remove the cap. If difficulty is experienced stop, let vein rest for a minute and try removal again very slowly. Gentle flushing of the catheter with saline may assist in guidewire removal.

Fig. 2 (Puncture with breakaway needle, peelable cannula or Neocath Split): Insert the catheter through the needle/cannula, using non-toothed forceps. Advance the catheter using short, steady strokes.

(Feeding the catheter with an anatomical forceps): If difficulty is experienced advancing the catheter, gentle flushing of the catheter with

saline solution may serve to stiffen the catheter slightly and may also assist in opening venous valves which may be obstructing advancement. When the final position is reached, stabilize the catheter by pressing lightly on the skin, 2-3 inches above the insertion site and withdraw the needle/cannula over the catheter.

INTRODUCER REMOVAL

The introducer is removed from the catheter in five stages: Stabilize the catheter to prevent withdrawing the catheter with the needle.

Important Caution: AT NO TIME SHOULD THE CATHETER BE WITHDRAWN BACK THROUGH A SPLITTING NEEDLE. IF IT BECOMES IMPOSSIBLE TO ADVANCE THE CATHETER INTO A SATISFACTORY POSITION, THEN THE NEEDLE AND CATHETER MUST BE WITHDRAWN SIMULTANEOUSLY. THE RESULT OF WITHDRAWING A CATHETER BACK THROUGH THE NEEDLE CAN BE CATHETER EMBOLISM.

Fig. 3 (withdraw the introducer from the vessel): Do not split or peel away the needle until it is completely removed from the vessel.

Fig. 4 gently pinch the wings of the needle completely together to break the bottom and the top part of the needle hub.

Fig. 5: Peel the wings of the needle apart, until the needle has peeled to within app. 1 mm of the tip. The top part of the needle will be completely open, and the catheter can be carefully removed with forceps. (Note: If considered preferable, the needle can be peeled completely into two pieces, by increasing, slightly, the peeling force, when resistance is felt, after peeling the needle to 1 mm from the tip.)

Anchor the catheter, using the fixation wings. If a vein in the arm is used, a small loop should be left in the catheter to prevent kinking of the catheter if the patient flexes or bends the arm.

IMPORTANT: ARRANGE FOR A CHECK X-RAY TO CONFIRM CORRECT PLACEMENT PRIOR TO STARTING ANY COURSE OF TREATMENT VIA THE CATHETER.

If the catheter is not to be used immediately, or its use is temporarily discontinued, it should be heparin locked to prevent clotting or blockage.

CAUTION: Do not exceed a bolus injection pressure of 1,5 bar. Do not use on an infusion pump that has infusion pressures that can exceed 1,5 bar as this will burst the catheter.

Check and record catheter position before starting an infusion, with the help of an x-ray photograph or during placement with the help of an intravascular electrocardiographic lead. The catheter tip must not be advanced into the heart (right atrium). The location of the catheter within the heart may cause cardiac tamponade, or cardiac arrhythmias.

To ensure utilisability of the catheter throughout the entire time period, it is necessary to make checks of the catheter position at regular intervals and to take care that the infusion pressure never exceeds 1bar. It is also necessary either to leave an infusion connected all the time, or to lock the catheter with heparin.

This device should only be used by, or under the instruction and responsibility of, an appropriately qualified person who is familiar with the product and the clinical procedure.

Complications associated with percutaneous catheters include vessel wall perforation, cardiac tamponade, haemorrhage, pleural and mediastinal injury, air embolism, catheter embolism, thoracic duct injury, septicæmia, bacteraemia, thrombosis, and accidental arterial puncture.

CAUTION:

Do not use this product for monitoring, diagnosis, control or treatment of a defect of the heart or the central venous system.

Do not bend the catheter or the extension line permanently to avoid damage to the catheter.

Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or bacterial infections.

Do not over stretch the catheter as it may rupture, and rebound into the insertion site, causing a catheter embolism.

Do not use the catheter for more than 29 days.

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