Laboratoire pharmaceutique VYGON

95440 Ecoeu - France

Tél. : +33 (0) 1 39 92 63 63 - Fax : +33 (0) 1 39 92 64 44

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SITIMPLANT / MINI-SITIMPLANT / MICRO-SITIMPLANT (for intravenous therapy)

Directions for use

PRODUCT DESCRIPTION

Vygon has developed an implantable port to overcome the problems associated with repeated venous access in patients needing long-term intravenous therapy. The Sitimplant is an entirely implantable port which is comprised of :

- A titanium reservoir. The area of silicone which covers the hollow of the reservoir is termed the septum and this is a self-sealing membrane which will be pierced to introduce medications, to feed intravenously and to take blood samples.
- a radio-opaque silicone (or polyurethane) catheter which is introduced either by surgical cut-down (cephalic or jugular veins) or by a percutaneous technique into the subclavian or jugular veins and then usually positioned so that its tip is in the superior vena cava. Depending upon the model of Sitimplant the catheter is either an integral part of the port (pre-attached) or the catheter is a separate component which is attached to the port during surgical placement.

N.B. Priming volume/Dead space : Sitimplant : 0.51 ml - Mini Sitimplant : 0.29 ml - Micro Sitimplant : 0,17 ml

•Silicone venous catheter with :

- *an internal diameter of 0.8 mm*

For a dead space of 60 cm : 0.30 ml - For a dead space of 10 cm : 0.05 ml

- *an internal diameter of 1.1 mm*

For a dead space of 60 cm : 0.60 ml - For a dead space of 10 cm : 0.10 ml

- *an internal diameter of 1.6 mm*

For a dead space of 60 cm : 1.10 ml - For a dead space of 10 cm : 0.20 ml


- *an internal diameter of 2.0 mm*

For a dead space of 60 cm : 1.88 ml - For a dead space of 10 cm : 0.31 ml


Example :

Mini-Sitimplant + 30 cm length of silicone catheter (with an i.d. of 1.1 mm)
 $0.29 \text{ ml} + (0.10 \text{ ml} \times 3) = 0.59 \text{ ml}$


Sitimplant / Mini-Sitimplant / Micro-Sitimplant in Cross-section


Connecting the catheter to the port

- Position the anti-kink sleeve on the proximal end of the catheter

- Attach the catheter to the connector on the port,
- taking care to ensure that the catheter is advanced fully over the barbs of the connector

- Bring the anti-kink collar over the port/catheter junction and make sure that the proximal side of the collar is fully against the wall of the port.

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PRIOR TO PLACING A SITIMPLANT

An implantable port requires surgical placement and the same procedures will apply as for any pre-operative patient.

- suitability of the patient for surgery determined by the patient's medical history physical examination and the completion of appropriate blood tests, X-rays etc.
- The area chosen for the implantation being shaved, washed and the skin surrounding the site being disinfected.
- Prepare the patient by explaining the procedure and the benefits and advantages to him of an implantable port.

N.B. Before handling the catheter all traces of talcum powder should be removed from the sterile gloves by washing them with sterile water. Use a Hubsite or an Injectsite to flush the catheter and port with a heparinised I.V. solution.

PLACEMENT OF THE PORT

The procedure is comparable to implanting a permanent pace-maker and it is usually carried out under local anaesthesia.

The choice of implantation site will depend upon the type of therapy and will be agreed with the patient by the physician, surgeon or anaesthetist.

There are two methods of implanting a Sitimplant:

The first of these is **surgical placement** which involves a "cut-down" into jugular or cephalic vein and the introduction of the catheter via the vein into the superior vena cava. The other end of the catheter is then tunnelled to the intended implantation site and is cut to the required length before being connected to the port. The port is then placed in a subcutaneous pouch that has been previously prepared in the chest wall.

The port and catheter are flushed with a heparin solution and the patency of the assembly confirmed before the port is sutured into place and the pouch closed.

The second method of implanting a port is termed "**percutaneous**" and uses a sheath-introducer (Desilet) to catheterise the vein.

Warning :

If the subclavian vein is used for catheterisation then care should be taken to avoid the catheter entering the subclavian vein medially because such placement can lead to compression of the catheter between the first rib and the clavicle and to severing or damage to the catheter. A fluoroscopic or radiographic confirmation of catheter placement should be made to ensure that the catheter is not being pinched between the first rib and the clavicle.

- If the catheter is an integral part of the port (pre-attached) this will necessitate the use of a "peel-apart" type of "Desilet" sheath-introducer. If the catheter is separate from the port then the introducer-sheath may be either of the "peel-apart" type (Vygon code 1146) or the standard type of sheath-introducer (Vygon code 1129 or 1137).

a) Percutaneous method using a standard sheath-introducer (Vygon code 1129 or 1137).

The Desilet is first placed in the subclavian or internal jugular vein and the catheter then introduced through the sheath of the Desilet. The sheath is then removed and the catheter tunneled into position before being cut-to-length and connected to the Sitimplant which is in the previously prepared subcutaneous pouch.

b) Placement using a "peel-apart" Desilet (Vygon code 1146).

The first step is to prepare the subcutaneous pouch. The Sitimplant with the pre-connected catheter is placed into the pouch. The catheter is tunneled, cut to the pre-determined length and then introduced into the Desilet which is positioned in the chosen vein. Peeling the sheath open while progressively advancing the catheter into the vessel allows the sheath to be removed completely with the catheter in place.

Regardless of the method of introducing the catheter, the correct positioning of its distal tip should be confirmed radiologically and, in the event of difficulty for whatever reason, diluted contrast media should be injected into the catheter to establish the tip position with certainty. **The distal tip of the catheter should be positioned in the superior vena cava and 2 cms above the right atrium.** The port and catheter should be carefully flushed and then heparinised. When suturing the subcutaneous pouch, assure that the suture is kept at a distance from the septum. After placement of the port sutures can normally be removed after 10 days. Unless the port has been used immediately after placement, the heparinisation carried out during placement will suffice for the 10 day period until the removal of the sutures.

MONITORING AND CARE OF THE SITIMPLANT

After each use of the port it will need to be flushed and then heparinised to prevent clot formation. This applies whether the port is accessed for an injection, an infusion, a blood sample or a blood transfusion. If the port is not being used then a monthly heparinisation is recommended.

Observe the following rules when using Sitimplant :

- Use only 10 ml (or larger) syringes
- Follow strict aseptic technique

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- Flush the port between each use and between each substance administered via the port
- Never use force to inject or infuse into the port.

ROUTINE PROCEDURES PRIOR TO USING THE PORT

Before accessing the port :

- Palpate the area over the port to identify and locate the port septum
- Check for any inflammation or reddening of the skin surrounding the port and for any sensitivity or pain in the area
- Carry out a preliminary disinfection of the area, at this stage it is not necessary to wear sterile gloves.

INJECTION VIA THE PORT

You will need :

- 3 10 ml syringes
- 1 3-way stopcock
- 1 Hubsite (Vygon code 1239) or an Inject-site
- 3 ampules of saline solution
- The prescribed medication
- Heparin solution (500 u/ml for adult - 100 u/ml for children)

Procedure

- WASH HANDS THOROUGHLY
- Put on sterile gloves
- Carry out a second disinfection of the area around the port
- Prepare the syringes : saline solution for flushes, heparin solution and medication
- Connect the syringe containing the saline solution to the 3-way stopcock and the Hubsite (or Inject-site)
- Prime the whole assembly and close the stopcock
- Anchor the port between the thumb and the first two fingers of the non-dominant hand.
- Insert the needle perpendicularly to the port septum taking care to avoid any lateral movement of the needle
- Advance the needle through skin and septum until contact is felt between the metal needle and the bottom of the reservoir
- Open the 3-way stopcock and check for venous blood return
- Flush with saline solution. No resistance should be felt to the flush and if there is any resistance then this indicates an obstruction of the catheter. In the event of an obstruction never attempt to clear it by forcibly injecting into the port as this

could rupture or split the catheter. During flushing, check that there is no swelling around the puncture site. If swelling is detected at the site when the port is flushed then remove the needle and stop the injection immediately

- Close the 3-way stopcock
- Connect the syringe(s) with the medication(s) and flush the port between each type of medication
- Flush again when all the injections are completed
- Finally, heparinise the Sitimplant with 3 ml of heparin solution (500 u/ml adult - 100 u/ml child)

CAUTION

If more than one drug is given during access to the port then it is essential to flush the port with saline solution after injection of each drug to minimise drug incompatibility reactions.

REMOVING THE NEEDLE

After flushing the port with saline solution the final step will be to heparinise the port and remove the needle. The two essential precautions when removing the needle are to stabilise the port by pressing it down with two fingers as force is applied to pull the needle out and to continue to inject the heparin solution whilst the needle is being withdrawn. Maintaining a positive pressure on the syringe plunger as the needle is withdrawn, prevents air being aspirated into the catheter and the attendant risk of clotting. After removal of the needle the site may be covered with a sterile, adhesive transparent dressing such as "Dermafilm".

N.B. The use of the accessory "Digiprotect®" (Vygon REF 701) will protect the user against needle-stick injuries during needle removal.

CONTINUOUS INFUSION

The method is identical to that for an injection except that a right-angled Hubsite needle is recommended for patient comfort and ease of dressing. A sterile dressing should be used and the entire site protected by a transparent, semi-permeable adhesive dressing such as "Dermafilm".

If the infusion lasts several days then the dressing should be renewed every two-three days.

All continuous infusions should be made via a delivery system that has the safety of a flow-control device.

BLOOD SAMPLING

Sampling via the Sitimplant uses the same equipment as injecting, only the flushing and the heparinisation sequence is different.

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- Flush the port with 5-10 ml of saline solution
 - Aspirate an initial 10 ml of blood which will normally be discarded (when aspirating, remember to do so progressively and steadily because a strong sudden negative pressure can cause the catheter walls to collapse)
 - Aspirate the amount of blood required for the tests
 - Flush with 20 ml of saline solution
 - Complete a final heparinisation (see notes)
- N.B. The 3-way stopcock must be closed to isolate the port whenever the syringes are changed.

POSSIBLE COMPLICATIONS

These can mostly be avoided by observing the recommendations for use.

No blood return but injection is possible : If it is possible to inject 20 ml of saline solution into the port without force and if this does not provoke swelling or pain at the site then the port can be used.

No blood return and injection is not possible : Check that the needle is correctly positioned in the port because failure to advance the needle tip completely through the septum will prevent both injection and aspiration.

Obstruction or blockage of the port or catheter as a result of clotting can normally be resolved by administering a fibrinolytic (see below). This is not the case when clotting is caused by the precipitation resulting from a reaction between incompatible substances having been administered through the port. If this occurs there is no alternative but to remove the port.

Skin/tissue necrosis is usually linked to an accidental injection of a chemotherapy agent which has been made subcutaneously instead of into the port. Use of the port should be stopped and the patient kept under close medical observation.

Infection levels are usually lower with an implanted port than with a conventional percutaneous central venous catheter. In addition to signs of sepsis around the port itself, infection should also be suspected in the event of unexplained feverishness or a temperature spike in a patient with an implanted port. Blood cultures should be made from the port area and from a peripheral vein. Normally an aggressive antibiotic treatment (systemic and also possibly localised) will overcome the infection and allow the port to be kept in place.

The risk factors for **vascular thrombosis** with an implanted port are no different from those with a conventional, percutaneous central venous catheter.

CAUTION :

Do not allow the use of syringes smaller than 10 ml. This is to avoid the high pressures generated by small syringes.

For information :

DISINFECTION : Povidone-iodine
or Chlorhexidine

HEPARINISATION : 3 ml of heparine solution at 500 U/ml for adults and
3 ml at 100 U/ml for children

DE-CLOTTING
WITH UROKINASE Use a 1000 U/ml solution
Inject 1 to 2 ml WITHOUT FORCE and wait
1 to 2 hours before re-aspirating the Urokinase, flushing and then heparinising

DRESSING Apply local pressure for 2 minutes and then apply a
"Dermafilm"
(code 38.15.10)

Pyrogen free. Sterility guaranteed unless unit pack has been damaged or opened.

Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or bacterial infections.

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