CATETERI OMBELICALI

Codici 270. / 1270. / 1272. / 1274.

Contenuto del kit:

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Code 270. / 1270. / 1272. / 1274.

Content of the kit:

270.xx: PVC single lumen umbilical catheter, 40cm long

1270.xx: Polyurethane single lumen umbilical catheter, 2,5Fr 30cm, other sizes 40cm long

1272.xx: Polyurethane double lumen umbilical catheter, 20cm long

1274.xx: Polyurethane double lumen umbilical catheter, 40cm long

Indications: The catheter is used short term catheterisation through the umbilical vessels.

Single lumen catheters can be used for catheterisation of the umbilical vein or arteries. The double lumen catheters can only be used for insertion into the umbilical vein.

Contra-indications: Umbilical catheters must not be used in case of peritonitis, necrotizing enterocolitis, omphalitis, omphalocele, evidence of local vascular compromise in lower limbs or buttocks area, or portal venous hypertension.**Insertion technique:** Use strict aseptic technique when handling and introducing the product. Disinfect carefully the umbilical cord and the abdominal wall. Place a sterile drape with central opening.

1. Elevate the cord vertically using toothed forceps. Cut the cord horizontally 1 cm above the skin with a scalpel blade. Remove the clots which may obstruct the vessel lumen. Control any bleeding by gentle tension on the umbilical tape.

2. Grasp the stump of the cord and hold it upright. Identify the vessels.

3. Prime the catheter with saline solution.

4. Carefully dilate the vessel in which the catheter will be inserted.

5. Insert the catheter into the vessel with forceps and advance the catheter directing it toward the pelvis.

6. Use the cm markings on the catheter to note the length of catheter inserted.

7. Place a suture in the wall of the cord taking care not to puncture the catheter.

8. Dress the catheter according to your local policies.

Control of positioning: Always check the location of the catheter by radiography.

For arterial catheterisation, the catheter tip should be placed in the lower abdominal aorta below the origin of the renal and inferior mesenteric arteries and above the bifurcation of the aorta (i.e., at the bottom of, or just below, the 3rd lumbar vertebra).

For venous catheterisation, the catheter tip should be placed beyond the ductus venosus in the central venous system (inferior vena cava)

Note: in emergency situations, the catheter can be inserted 2-3 cm (until blood is returned) and emergency meds can be given. The catheter should not be left in this placement and should be removed after giving emergency medications.

Nursing care: The catheter and insertion site should be inspected daily for any indications of a developing complication. Signs and symptoms at the insertion site include warmth, redness, swelling, pain or tenderness, circulatory compromise, discharges... If any of these are present, notify the appropriate doctor and record in the nursing records. The catheter's general condition should also be noted.

Flush the catheter before and after each use with a saline solution.

Catheter Cleaning/Dressing: Catheter dressing policy should follow the protocol adopted by the institution involved.**CAUTION:** Be aware that organic solvents such as alcohol or acetone may interact with PUR catheter material and weaken it.**Removal:** The catheter insertion site should be cleaned carefully prior to catheter withdrawal.

The catheter is removed by gentle sustained traction close to the exit site (2-3 cm). Do not over stretch the catheter.

After removal, inspect the catheter to confirm its integrity.

Potential complications: Complications associated with umbilical catheters include ischemia, thrombosis, vasospasm, haemorrhage, infection, vascular perforation, air embolism...

This device should only be used by, or under the instruction and responsibility of, an appropriately qualified person who is familiar with the product and the clinical procedure.

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