EPICUTANEO-CAVE-CATHETER/ EPICUTANEO STAR / EPICUTA- NEO-CAVE-CATHETER EXPERT

Instructions for Use

Art.-No. 1284., 2184., 6184., 8184.

Important cautions:

At no time totally unscrew the blue hub.

At no time withdraw the catheter back through the winged needle.

At no time exceed the working or bolus pressure limits of 1 bar/760mm Hg and 1,2 bar/900mm Hg. N.B.: Art.-No. 1284 catheter is manufactured using 24 gauge XRO polyurethane tubing. Art.-No. 2184, 6184 and 8184 catheters are manufactured using 24 gauge XRO silicone tubing. This is a delicate material which should be treated with care at all times. The catheter has a detachable self-compression hub which grips onto the metal tube which is already in position inside the proximal end of the catheter.


Introduction:

Selecting a vein

This catheter is suitable for insertion into peripheral veins. Recommended insertion sites are median cubital veins at the elbow and wrist, or the basilic vein. The long saphenous veins of the ankle may also be used. This catheter is not suitable for insertion through non-superficial veins.

Insertion technique


1. The introduction of this catheter should be performed using strict aseptic technique. The site of needle insertion should be cleaned thoroughly using an appropriate prep solution. Before handling the product we recommend that latex gloves should be washed to remove all residue of starch or glove powder (Fig. A).


2. Determine the length of catheter needed to reach the final position, using a measuring tape. Allow for an additional loop of catheter for external fixation. If required shorten to length by cutting only from the distal end to prevent accidental removal of the metal tube. Be careful to make a clean, vertical cut. The length of a shortened catheter should be recorded prominently in the patient's notes.

3. Prior to checking catheter patency tighten compression hub. Check the patency of the catheter and extension tube with a fluid-filled syringe. Do not occlude the catheter during this procedure.

4. Using the 19G winged needle, perform venipuncture. Use only the winged needle provided (Fig. B).

5. Carefully insert the catheter, distal tip first, through the winged needle. Use only nontoothed forceps (Fig. C).

6. Partially unscrew the compression hub and remove the catheter. Do not totally remove the blue portion (Fig. D+E).

7. Advance the catheter carefully through the needle to the final position. Do not at any time withdraw the catheter back through the winged needle. This may cause puncture or embolisation of a portion of the catheter.


8. Carefully withdraw the winged needle, taking care to avoid pulling on the catheter. Fix the catheter in the vein during needle removal by maintaining gentle pressure over the vein distal to the needle bevel (Fig. F).

Take particular care passing the needle over the metal insert on the proximal portion of the catheter (Fig. G).

Partially unscrew the compression hub and remove the catheter. Do not totally remove the blue portion.

9. Insert the proximal end of the catheter containing the metal tube into the blue portion of the compression hub until a definite stop is experienced. Tighten the blue portion. The thick black ring on the catheter should not be visible at the distal end of the blue compression hub (Fig. H+I).

10. Check the connection between the compression hub and the catheter by pulling gently.

Caution: Do not use small syringes as these can generate very high pressures. It is possible to generate 4 or 5 times the maximum safety pressure, with any size of hand held syringe. Subjecting the catheter to pressure above 1,2 bar can result in catheter rupture and embolism.

11. Flush catheter at this point without exceeding pressure limits (see Important Cautions). Do not attempt to aspirate blood to determine location.

12. Secure the catheter in a loop without causing any strain of the tubing. Do not apply adhesive tape/dressings directly onto the tubing as this can damage the delicate tubing when removed. The blue hub and extension line should be attached to the patient to prevent any strain on the catheter.

13. Obtain and view a check x-ray to establish correct positioning prior to commencing an infusion of hypertonic solution.

Nursing care:

The catheter and insertion site should be inspected daily for any indications of a developing complication. Signs and symptoms at the insertion site include warmth, redness, swelling, pain or tenderness and discharges. If any of these are present, notify the appropriate doctor and record in the nursing records. The catheter's general condition should also be noted.

Catheter Cleaning/Dressing

Catheter dressing policy should follow the protocol adopted by the institution involved.

CAUTION: Be aware that organic solvents such as alcohol or acetone may interact with PUR catheter material and weaken it. Alcohols such as ethanol and isopropanol do not change the physical properties of the silicone elastomer if they are used briefly without prolonged contact and if the silicone catheter is then dried. The use of organic solvents such as acetone, chlorinated solvents, ethyl acetate, etc. is inadvisable. Do not expose this catheter to contact with iodine containing compounds! These products can degrade or irreversibly deform the silicone elastomer.

Removal

The catheter insertion site should be cleaned carefully to remove all sticky residue prior to catheter withdrawal.

The catheter is removed by gentle sustained traction close to the exit site (2-3 cm). Do not over stretch the catheter as it may rupture and rebound into the vein causing a catheter embolus. A finger should be placed directly over the vein (without applying pressure) so that the catheter can be fixed immediately should catheter rupture occur.

If the catheter remains firmly attached despite sustained gentle traction then surgical removal is needed.

Important cautions in use

1. Do not use infusion equipment which can exceed the working pressure of 1,0bar max/760 mm Hg (see maximum flow capability tables). Only use infusion equipment complying with standards, which do not exceed a shut-off pressure of 1,0 bar.

2. Bolus injections should be slow and must not exceed the maximum bolus pressure of 1,2bar/900 mm Hg.

3. Do not use syringes smaller than 10 cc. Smaller syringes generate higher pressures than larger ones.

4. Check compression hub regularly to ensure black ring is not visible (see note 9) and blue hub for tightness.

5. A running infusion should be maintained at all times to ensure patency. A heparin lock may be substituted for an infusion.

6. If a heparin lock is used, ensure that the catheter, extension tube and hub assemblies are flushed to exclude all drugs, TPN solutions, air, etc. from the system prior to instilling the heparin lock.

7. The extension tube must not be clamped after heparin locking and then released. This will cause blood to be aspirated back into the catheter lumen. This can result in a blocked catheter.

8. The clamp provided on the extension tube is for temporary clamping only during connections and disconnections. Prolonged clamping of the extension tube will cause damage to the tubing.

Complications associated with percutaneous catheters include vessel wall perforation, cardiac tamponade, haemorrhage, pleural and mediastinal injury, air embolism, catheter embolism, thoracic duct injury, septicaemia, bacteraemia, thrombosis, and accidental arterial puncture.

Do not use this product for monitoring, diagnosis, control or treatment of a defect of the heart or the central venous system. Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or bacterial infections.

Do not bend the catheter or the extension line permanently to avoid damage to the catheter.

Do not use the catheter for more than 29 days.

This device should only be used by, or under the instruction and responsibility of, an appropriately qualified person who is familiar with the product and the clinical procedure.

Caution: Federal U.S. law restricts this device to sale by or on the order of a physician.

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