



Value Life

IV Accessories
Solutions for safety



Blunt fill needles

Including blunt fill needles
with 5 micron filter


aim at **ZERO**

Now in
Tyvek
packaging

Help reduce the risk

of needle-stick injuries to healthcare workers and patients.

Compared with sharp-bevelled hypodermic needles, blunt fill and blunt filter needles help **reduce the risk of needlestick injuries**¹ to healthcare workers when drawing-up and preparing medication in compliance with EU Directive 2010/32/EU.



Colour-coded hub for easy identification

- Blunt fill needle (BFN) for use with vials.
- Blunt fill needle with 5 micron filter (BFNF) for use with glass ampoules.

Ease of use

- Both needles have a conventional luer needle hub, compatible with all standard syringes.

Reduction in medication contamination

- The 5 micron filter in the hub of the BFNF means any glass particles larger than 5 microns are filtered out preventing them from being drawn into the syringe prior to injection.

The use of a 5 micron filter needle or straw to withdraw medication from the ampoule can reduce the number of particles aspirated. ^(2,3)

How to use...

Step-by-step guidelines for our blunt fill needles.

Step 1

- Peel open the packaging and connect the syringe.



Step 2

- Remove the needle guard.

Step 3

BFN (single-use only)

- Insert the needle into the centre of the vial stopper at a 90° angle. This will help reduce the risk of rubber fragments from the side of the vial stopper contaminating the medication.



BFNF

- Draw-up the medication from the glass ampoule.

Step 4

- Once the medication has been drawn-up and is in the syringe ready for administration, remove the BFN / BFNF and dispose of in a sharps container.



Glass particle contamination can occur when opening single-use glass ampoules of medication, and injection of these particles has been associated with phlebitis, pulmonary thrombi or microemboli, and end organ granulomas or inflammation. ^(2,3)

Product Specification

Vygon Code	NHSSC Code	Description	Quantity
110022	FTR1922	18G x 1½"	100
110022F	FTR1923	18G x 1½" with 5 micron filter	100

Supplied in DuPont™ Tyvek® packaging

The unique structure of our packaging creates a tortuous path for a superior microbial barrier and excellent strength properties. Made of high-density polyethylene (HDPE), Tyvek® offers all the best characteristics of paper, film and fabric in material. This unique balance of properties makes our packaging lightweight yet strong; vapour-permeable, yet moisture and chemical resistant, as well as puncture and tear resistant.

References

1. Adams D, Elliott T. Impact of safety devices on occupationally acquired needlestick injuries: a four-year prospective study. *Journal of Hospital Infection* 2006; 64: 50–5
2. Preston ST, Hegadoren K. Glass contamination in parenterally administered medication. *J Adv Nurs*. 2004;48(3):266-70
3. Zahir AF, Choy CY, Rushdan R. Glass particle contamination of parenteral preparations of intravenous drugs in anaesthetic practice. *Southern African Journal of Anaesthesia and Analgesia*. 2008;14(3):17-9.

All references are available on request.

For further information, please contact: vygon@vygon.co.uk

The specifications shown in this leaflet are for information only and are not, under any circumstances, of a contractual nature. This document is intended for use in the UK only.

Vygon (UK) Ltd • The Pierre Simonet Building • V Park • Gateway North
• Latham Road • Swindon • Wiltshire • SN25 4DL

Tel: 01793 748800 • Fax: 01793 748899 • Twitter: @vygonuk

Web: www.vygon.co.uk Content correct as of: 06/2017

Code: DXJB0100015 v1 Copyright Vygon (UK) Ltd 2017



www.vygon.co.uk