

# INTRAVASCULAR THERAPIES

ULTRASOUND

**vysion** Andy

PRODUCT GUIDE



Value Life

## WHY CHOOSE VYGON?




**OVER  
60  
YEARS  
EXPERIENCE**

## **PATIENT CARE** LIES AT THE HEART OF WHAT WE DO



**CLINICAL  
EXPERTISE**



WE ARE A  
**CARBON  
NEUTRAL  
COMPANY**

MORE THAN  
**100 MILLION**  
PRODUCTS IN  
**110 COUNTRIES**



## CONTENTS

Pocket ultrasound system for vascular access	3
Ergonomical design	4
Quick and easy access	4
Other advantages	5
Midlines	6
Our commitment to the environment	7
Product range and ordering codes	7

## OUR SERVICE OFFERING



### **Customer Service & Technical Support**

Talk through your enquiry with our dedicated Teams.



### **Education & Training**

We offer our customers a variety of valuable and comprehensive training options to help you and your teams meet your training requirements.

# HANDHELD ULTRASOUND SYSTEM FOR VASCULAR ACCESS

## ULTRASOUND SYSTEM WITHIN ANY REACH

Looking for an handheld or mobile ultrasound system offering:

- High quality imagery
- Easy mobility
- Compatible with any device: PC, tablet and even smartphone!

vysion Andy can assist you with all your vascular access procedures from echo-tracking and ultrasound-guided puncture, to tip location.



### Small Linear probe



Small: 170mm\*56mm\*22mm



Super light: only 220g




Central frequency: 7.5MHz  
Frequency range: 6.0-11MHz

### Other applications:

Nerve, Anaesthesia, Vascular, Small Parts, MSK, Orthopaedics Available in Black & White or in Colour mode.

# ERGONOMIC DESIGN

Advanced workflow thanks to an easy to operate probe



# QUICK AND EASY ACCESS TO THE VEIN FROM ANY DEVICE

Vysion Andy enable you to perform scan exams anywhere using the application sonoiQ.



The probe has its own wifi signal to work, no need for an external one!



The wireless probe is compatible on



## OTHER ADVANTAGES:



**Waterproof** wireless probe convenient for sterilisation.



**Low energy consumption** (max. 3W) allowing 8 hours of continuous scanning time.



## VYSION ANDY: MAKE CATHETER INSERTION EASY!

Vysion Andy wireless probe allows you to have access to ultrasound for any type of catheter you can use for vascular access.

Midline




Port



PICC

# smartmidline

Simplicity of insertion



With the increased demand for ultrasound-guided placement, vysion Andy is the tool that allows you to place midlines better by visualising the vein and structure of the arm.

Place more in the upper arm and avoid the complications of forearms' smaller veins, lower dwell time and possible dislodgement with our portable handheld ultrasound device.

## OUR COMMITMENT TO THE ENVIRONMENT

Vygon UK has been working towards creating a sustainable future for many years and we are a Carbon Neutral Company. We believe that sustainability isn't just about meeting current needs, but more importantly, it is about ensuring we are here for the long term and are paving the way for a bright tomorrow. We strive to deliver enhanced Corporate Social Responsibility (CSR), ensuring we manage the social, economic, and environmental effects of Vygon UK's operations responsibly in line with public expectations.



**DISCOVER OUR JOURNEY**  
scan with your smart device



## ORDERING CODES

Vysion Andy	
Product code	vysion Andy version
VNY1YO8H03301UK	VYSION ANDY WIRELESS LINEAR PROBE COLOR + DICOM
VNY1YO8H04302UK	VYSION ANDY WIRELESS LINEAR PROBE BW + DICOM



12/2023 - DXJ00000000 V01644 v1

FOR FURTHER INFORMATION, PLEASE CONTACT: [info@vygon.co.uk](mailto:info@vygon.co.uk)

The specifications shown in this leaflet are for information only and are not, under any circumstances, of a contractual nature. This brochure has been printed on responsibly sourced and sustainable material. To help us reduce our carbon footprint all of our literature is available electronically either from your Product Specialist, on our website or by emailing [info@vygon.co.uk](mailto:info@vygon.co.uk)

VYGON (UK) LTD, THE PIERRE SIMONET BUILDING, V PARK,  
GATEWAY NORTH, LATHAM ROAD, SWINDON, WILTSHIRE SN25 4DL  
RECEPTION: +44 (0)1793 748800 [WWW.VYGON.CO.UK](http://WWW.VYGON.CO.UK)

 [@vygonuk](https://www.instagram.com/vygonuk)

